

## M2培养基

产品货号: T10512

产品规格: 100ml/500ml

保存: 2-8°C保存, 有效期1年。

### Product Description:

M16 and M2 are modified Krebs-Ringer bicarbonate solutions, which are similar to Whitten's Medium. M16 Medium is one of the most common media for in vitro culture of preimplantation stage embryos. M16 contains pyruvate and lactate as energy sources since preimplantation embryos do not utilize glucose efficiently. M2 Medium is a further modification of M16 that substitutes HEPES buffer in place of some of the bicarbonate. M2 is used for collecting and handling embryos for prolonged periods outside a CO<sub>2</sub> incubator.

### Formulation: Ready-to-use, sterile filtered solution:

Component	Concentration (g/L)
CaCl <sub>2</sub> •2H <sub>2</sub> O	0.251
MgSO <sub>4</sub> (anhydrous)	0.165
KCL	0.356
KH <sub>2</sub> PO <sub>4</sub>	0.162
NaCl	5.532
NaHCO <sub>3</sub>	0.35
Albumin, Bovine Fraction V	4.0
D-Glucose	1.0
HEPES	5.43
Phenol Red•Na	0.01
Pyruvic Acid•Na	0.036
Lactic Acid•Na	4.35
Specifications	
pH	7.4±0.3
Osmolality (mOsm)	285±25

### References:

1. Whitten, W.K. (1971) Embryo medium: Nutrient requirements for the culture of preimplantation embryos in vitro. *Adv. Biosci.*, 6: 129-141.
2. Quinn, P., Barros, C., and Whittingham, D.G. (1982) Preservation of hamster oocytes to assay the fertilizing capacity of human spermatozoa. *J. of Reprod. Fertil.*, 66: 161-168.
3. Hogan, B., Costantini, F., and Lacy, E. (1986) *Manipulating the Mouse Embryo*. Cold Spring Harbor Laboratory. Cold Spring Harbor, NY. Section G, pp. 249-252.



扫一扫 加微信

郑州乐业生物科技有限公司  
Zhengzhou Leye-Bio Biotechnology Co.,Ltd

地址:郑州市高新区红松路36号龙鼎企业中心一期1号楼5楼25号

免费电话: 400-611-0007 13671551480 13643719799

Q Q: 807961520 731791866

邮箱: zzlybio@126.com